Directions for use



Ambu[®] **Redi-ACE Adult**

Ambu[®] **Redi-ACE Mini**

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
www.ambu.com

492-4800-40 - V01 - 03/2011 - Printed in China

Sizes	
1	Infant
2	Paediatric
3	Neckless
4	Short
5	Regular
6	Tall

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Warnings

The Ambu Redi-ACE collar is intended for application by licensed health care practitioners as directed by a physician or other medical authority. Permanent injury may result from use other than as directed by a physician or other medical authority.

An extrication collar is only part of total patient immobilization.

It is important that practitioners are trained and practiced in the use and application of the device and that all instructions are read and understood prior to application of the collar.

All directions are to be considered as guidelines and not an attempt to define medical practice. Use the device only after review and approval from the supervising medical authority.

The collar is intended for Single Patient use only due to the risk of cross contamination. It should be discarded after use in accordance to local regulations.

Check the device before use.

Patient to be monitored while device is in use.

1. Intended use

The Ambu® Redi-ACE collar is a one-piece rigid cervical spine immobilization device designed to assist the rescuer with the maintenance of neutral alignment, prevention of lateral (side-to-side) sway and anterior-posterior (forward and backward) flexion and extension of the cervical spine during transport and routine patient care or movement.

2. Specifications

Storage Temperature: -30°C ~ +70°C

Operational Temperature: -10°C ~ +40°C

3. Parts/Materials

Main panel: Polypropylene
Sliding panel: Polypropylene
Buttons: Polyamide

Velcro: Polyamide
Foam: IXPE

4. Sizing and Application

4.1 Sizing (1) (2)

- Measure the distance between an imaginary plane drawn horizontally and immediately below the patient's chin and second horizontal plane drawn immediately on top of the patient's shoulder.
- Compare this distance with the distance from the collar sizing line to the lower aspect of the plastic collar body (not the foam).
- Adjust the collar to the appropriate size.
- Simply pull the collar apart until the distance between the sizing line and the plastic collar body equals your finger measurement.
- Engage the safety locks by pushing DOWN the locking buttons.
- If the Ambu Redi-ACE needs to be resized, disengage the safety locks by lifting the locking buttons.
- The Ambu Redi-ACE may now be adjusted to the appropriate size.
- Now engage the safety locks by pushing DOWN the locking buttons.

4.2 Application (3)

- Proper application of any extrication collar requires two individuals. The first rescuer should maintain stabilization and neutral alignment of the head and neck.
- The second rescuer should slide the front of the collar along the patient's chest and position the chin piece. The collar body should rest on the top of the patient's shoulder and against the sternum without gaps. The patient's chin must rest securely in the chin piece in neutral alignment position.
- While holding the front of the collar in place, wrap the back of the collar around the back of the patient's head and neck and secure to the front of the collar with the Velcro® tab.
- Adjust the Velcro tab as necessary by holding the front cutout of the collar and tightening the Velcro tab.

4.3 Supine Application (4)

If the patient is supine, slide the back portion of the collar directly behind the patient's neck.

The collar may be applied from either side. Once the back portion of the collar is evenly positioned behind the head, the front portion of the collar can be positioned under the chin. While maintaining the proper positioning of the collar with one hand, attach the Velcro with the other hand. Ensure that the loop Velcro tab is firmly attached and parallel to the fixed hook Velcro tab on the front of the collar.

Warnings

Maintain neutral head and neck alignment during the entire procedure. Once the collar is properly applied and the patient is immobilized, the head should not be released. Further immobilize the head and neck with head blocks or a head immobilizer and secure the patient to the backboard with appropriate retention straps.

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